



Optional Refrigerated Storage Expiration Record

Refrigerated Storage Starting Date*	
--	--

***ATTENTION: GRAFT EXPIRES TWELVE (12) MONTHS FROM DATE RECORDED.**

Lipoderma™ is preserved frozen and **must be maintained at a temperature of -40°C or colder** excluding the alternate refrigerated storage condition listed below:

- **Refrigerated storage between 0°C to 10°C is acceptable for a maximum of 12 months from refrigeration start date.** Once placed at refrigerated temperatures, Lipoderma™ may NOT be refrozen, and the allograft must be used within twelve (12) months or discarded.
- If refrigerated storage is used, the facility must record refrigeration start date and establish a process to ensure the allograft is used within 12 months or discarded.

Ensure this record is maintained with the allograft.

Note: For allografts that reflect a (6) six-month expiration for refrigerated conditions on the package label, the 12-month expiration time frame applies.